


Anti-Nfkbib Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15653

Item Number: CM0027132



Introduction

Anti-NFKBIB Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Targets I-kappa-B-beta, a 38 kDa NF-kappa-B inhibitor that regulates cytoplasmic retention and nuclear translocation of NF-kappa-B.

[Learn More](#)

Parameter	Value
Product Name	NFKBIB Rabbit pAb
Remarks/Alias	IKBB; TRIP9; NFKBIB
Species	Human
Gene ID (Human)	4793
Gene ID	4793
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-356 of human NFKBIB (NP_002494.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15653
Protein Weight	38kDa
Shipping	Ice bag