

Anti-Nfkbiz Polyclonal Antibody For Wb And Elisa - Q9Byh8

Item Number: CM0014126



Introduction

Human NFKBIZ (I-kappa-B-zeta, IκB-zeta) rabbit polyclonal antibody validated for WB and ELISA, with cross-reactivity in rat. Ideal for studying NF-kappa-B regulation and inflammatory response pathways.

[Learn More](#)

Parameter	Value
Product Name	NFKBIZ Rabbit pAb
Remarks/Alias	IKBZ; INAP; MAIL; ikB-zeta; ikappaBzeta; IkappaB-zeta; I-kappa-B-zeta; NFKBIZ
Species	Human
Gene ID (Human)	64332
Gene ID	64332
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human NFKBIZ (NP_113607.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q9BYH8
Protein Weight	78kDa
Shipping	Ice bag