

Anti-Nfrkb Rabbit Monoclonal Antibody For Wb, Elisa - Q6P4R8

Item Number: CM0006500



Introduction

NFRKB Rabbit mAb for WB and ELISA detects human nuclear factor related to kappa-B-binding protein, a 139 kDa INO80 complex subunit involved in chromatin remodeling and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	NFRKB Rabbit mAb
Remarks/Alias	INO80G; NFRKB
Species	Human
Gene ID (Human)	4798
Gene ID	4798
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 370-495 of human NFRKB (NP_001137307.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6P4R8
Protein Weight	139kDa
Shipping	Ice bag