


## Anti-Nfu1 Polyclonal Antibody For Wb And Elisa - Q9Ums0

Item Number: CM0028496



### Introduction

High-quality rabbit polyclonal antibody against NFU1, suitable for Western blot and ELISA applications. Validated in human, mouse, and rat samples. Ideal for mitochondrial iron-sulfur cluster research.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | NFU1 Rabbit pAb                                                                                                     |
| Remarks/Alias    | Nfu; NifU; HIRIP; MMDFS; MMDS1; NIFUC; CGI-33; HIRIP5; NFU1                                                         |
| Species          | Human                                                                                                               |
| Gene ID (Human)  | 27247                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-230 of human NFU1 (NP_001002755.1). |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, ELISA                                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | Q9UMS0                                                                                                              |
| Protein Weight   | 28kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |