

Anti-Nfyb Rabbit Polyclonal Antibody For Wb, Elisa - P25208

Item Number: CM0027930



Introduction

NFYB rabbit polyclonal antibody for WB, ELISA. Reacts with human and mouse. Targets nuclear transcription factor Y subunit beta, a CCAAT box-binding protein involved in gene transcription regulation.

[Learn More](#)

Parameter	Value
Product Name	NFYB Rabbit pAb
Remarks/Alias	HAP3; CBF-A; CBF-B; NF-YB; NFYB
Species	Human
GeneID (Human)	4801
GeneID	4801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-207 of human NFYB (NP_006157.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P25208
Protein Weight	23kDa
Shipping	Ice bag