

Anti-Nfyc Rabbit Polyclonal Antibody For Wb And Elisa - Q13952

Item Number: CM0028655



Introduction

High-quality rabbit polyclonal antibody against NFYC (NF-YC, H1TF2A) for Western blot and ELISA. Validated in human and cross-reacts with rat. Ideal for transcription factor research and CCAAT-box binding studies.

[Learn More](#)

Parameter	Value
Product Name	NFYC Rabbit pAb
Remarks/Alias	HSM; CBFC; HAP5; CBF-C; NF-YC; H1TF2A; NFYC
Species	Human
GeneID (Human)	4802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-120 of human NFYC (NP_055038.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q13952
Protein Weight	50kDa
Shipping	Ice bag