

Anti-Ng2 Monoclonal Antibody For Flow Cytometry - Q6Uvk1

Item Number: CM0008892



Introduction

-conjugated rabbit monoclonal antibody against human NG2 (CSPG4) for flow cytometry. Validated for human samples. Useful for studying melanoma, angiogenesis, and cell migration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human NG2 mAb
Remarks/Alias	NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG
Species	Human
GeneID (Human)	1464
GeneID	1464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1583-2224 of human NG2 (NP_001888.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q6UVK1
Protein Weight	251kDa
Shipping	Ice bag