

Anti-Ng2/Cspg4 Rabbit Monoclonal Antibody For Wb, Elisa - Q6Uvk1

Item Number: CM0007015



Introduction

High-quality anti-NG2/CSPG4 rabbit monoclonal antibody validated for Western blot and ELISA. Targets human chondroitin sulfate proteoglycan 4 (251kDa), a key player in cell migration, melanoma invasion, and angiogenesis.

[Learn More](#)

Parameter	Value
Product Name	NG2/CSPG4 Rabbit mAb
Remarks/Alias	NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG; NG2/CSPG4
Species	Human
GeneID (Human)	1464
GeneID	1464
Immunogen	Recombinant protein (recombinant fusion protein containing a sequence corresponding to amino acids 1340-1563 of human NG2/CSPG4, NP_001888.2)
Source	Rabbit
Category	Monoclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q6UVK1
Protein Weight	251kDa
Shipping	Ice bag