

Anti-Ng2/Cspg4 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Fc, Elisa - Q6Uvk1

Item Number: CM0007381



Introduction

Recombinant rabbit monoclonal antibody targeting human NG2/CSPG4, suitable for WB, IF-P, IHC-P, FC, ELISA. Ideal for investigating cell migration, angiogenesis, and melanoma mechanisms; cross-reacts with mouse and human.

[Learn More](#)

Parameter	Value
Product Name	NG2/CSPG4 Rabbit mAb
Remarks/Alias	NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG; NG2/CSPG4
Species	Human
GeneID	1464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1583-2224 of human NG2 (NP_001888.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, FC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6UVK1
Protein Weight	251kDa
Shipping	Ice bag