

Anti-NG2/Cspg4 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q6Uvk1

Item Number: CM0022404



Introduction

High-quality anti-NG2/CSPG4 rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects human and mouse chondroitin sulfate proteoglycan 4, a key cell surface receptor involved in cell migration, angiogenesis, and melanoma progression. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	NG2/CSPG4 Rabbit pAb
Remarks/Alias	NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG; NG2/CSPG4
Species	Human
Gene ID (Human)	1464
Gene ID	1464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1950-2210 of human NG2/CSPG4 (NP_001888.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6UVK1
Protein Weight	251kDa
Shipping	Ice bag