

Anti-Ng2/Cspg4 Rabbit Polyclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Elisa - Q6Uvk1

Item Number: CM0027710



Introduction

KO validated rabbit polyclonal antibody against human NG2/CSPG4 (SWISS Q6UVK1). Suitable for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for cell migration, angiogenesis, and melanoma research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] NG2/CSPG4 Rabbit pAb
Remarks/Alias	NG2; MCSP; MCSPG; MSK16; CSPG4A; HMW-MAA; MEL-CSPG; G4
Species	Human
GeneID (human)	1464
GeneID	1464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1340-1563 of human NG2/CSPG4 (NP_001888.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UVK1
Protein Weight	251kDa
Shipping	Ice bag