


Anti-Ngb (Neuroglobin) Rabbit Polyclonal Antibody For Wb And Elisa - Q9Npg2

Item Number: CM0028939



Introduction

Rabbit polyclonal antibody against human neuroglobin (NGB), validated for WB and ELISA with cross-reactivity to mouse. Suitable for hypoxia, neuroprotection, and brain tissue research.

[Learn More](#)

Parameter	Value
Product Name	NGB Rabbit pAb
Remarks/Alias	NGB
Species	Human
Gene ID (Human)	58157
Gene ID	58157
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-151 of human NGB (NP_067080.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9NPG2
Protein Weight	17kDa
Shipping	Ice bag