

Anti-Ngdn Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nej9

Item Number: CM0013170



Introduction

Rabbit polyclonal antibody targeting human NGDN (Neuroguidin), recommended for WB and ELISA. Cross-reacts with rat. Suitable for studies on ribosome biogenesis, translation regulation, and neuronal development.

[Learn More](#)

Parameter	Value
Product Name	NGDN Rabbit pAb
Remarks/Alias	NGD; LCP5; CANu1; lpd-2; C14orf120; NGDN
Species	Human
Gene ID	25983
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-209 of human NGDN (NP_001036100.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8NEJ9
Protein Weight	36kDa
Shipping	Ice bag