


Anti-Ngdn Rabbit Polyclonal Antibody For Wb, Elisa - Q8Nej9

Item Number: CM0014704



Introduction

NGDN Rabbit pAb is a polyclonal antibody validated for WB and ELISA, with cross-reactivity to Human, Mouse, Rat. Ideal for studying ribosome biogenesis and translation regulation.

[Learn More](#)

Parameter	Value
Product Name	NGDN Rabbit pAb
Remarks/Alias	NGD; LCP5; CANu1; lpd-2; C14orf120; NGDN
Species	Human
GeneID (Human)	25983
GeneID	25983
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 61-209 of human NGDN (NP_001036100.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NEJ9
Protein Weight	36kDa
Shipping	Ice bag