

Anti-Ngf Polyclonal Antibody For Wb, Elisa - P01138

Item Number: CM0017084



Introduction

Rabbit polyclonal antibody against human NGF (Beta-nerve growth factor), validated for WB and ELISA, reacts with human and mouse. Suitable for studies of neuronal development and neuropathy. - P01138

[Learn More](#)

Parameter	Value
Product Name	NGF Rabbit pAb
Remarks/Alias	NGFB; HSAN5; Beta-NGF; NGF
Species	Human
GeneID (human)	4803
GeneID	4803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-241 of human NGF (NP_002497.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P01138
Protein Weight	27kDa
Shipping	Ice bag