


Anti-NGf Polyclonal Antibody For Wb, Elisa - P01138

Item Number: CM0015794



Introduction

Rabbit polyclonal antibody targeting human beta-nerve growth factor (NGF). Suitable for Western blot (WB) and ELISA; cross-reacts with mouse and rat. Ideal for neurobiology and growth factor signaling studies.

[Learn More](#)

Parameter	Value
Product Name	NGF Rabbit pAb
Remarks/Alias	NGFB; HSAN5; Beta-NGF; NGF
Species	Human
Gene ID (Human)	4803
Gene ID	4803
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 19-241 of human NGF (NP_002497.2))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P01138
Protein Weight	27kDa
Shipping	Ice bag