

Anti-Ngfrap1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q00994

Item Number: CM0029030

Introduction

NGFRAP1 Rabbit pAb validated for WB, IF/ICC, ELISA in human and mouse. Targets BEX3 (p75NTR-associated cell death executor) involved in apoptosis; ideal for neuronal death and zinc signaling research.

[Learn More](#)

Parameter	Value
Product Name	NGFRAP1 Rabbit pAb
Remarks/Alias	Bex; NADE; HGR74; NGFRAP1; DXS6984E
Species	Human
Gene ID (Human)	27018
Gene ID	27018
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human NGFRAP1 (NP_996798.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q00994
Protein Weight	13kDa
Shipping	Ice bag