


Anti-Ngly1 Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q96lv0

Item Number: CM0007437



Introduction

Rabbit monoclonal antibody against NGLY1 (Peptide-N(4)-(N-acetyl-beta-glucosaminy)asparagine amidase), validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat. Ideal for ERAD and glycoprotein degradation studies.

[Learn More](#)

Parameter	Value
Product Name	NGLY1 Rabbit mAb
Remarks/Alias	CDDG; PNG1; CDG1V; PNG-1; PNGase; NGLY1
Species	Human
GeneID (Human)	55768
GeneID	55768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NGLY1 (NP_060767.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96lv0
Protein Weight	74kDa
Shipping	Ice bag