


Anti-Ngly1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q96lv0

Item Number: CM0017164



Introduction

Rabbit polyclonal antibody targeting human NGLY1 (N-glycanase 1), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. Suitable for ER-associated degradation and glycoprotein quality control research.

[Learn More](#)

Parameter	Value
Product Name	NGLY1 Rabbit pAb
Remarks/Alias	CDDG; PNG1; CDG1V; PNG-1; PNGase; NGLY1
Species	Human
GeneID (Human)	55768
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NGLY1 (NP_060767.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96lv0
Protein Weight	74 kDa
Shipping	Ice bag