

Anti-Ngr3 (Rtn4Rl1) Rabbit Monoclonal Antibody For Wb, Ihc, If - Q86Un2

Item Number: CM0002550

Introduction

Recombinant rabbit monoclonal antibody targeting human Ngr3 (RTN4RL1), validated for WB, IHC, IF, and ELISA. Cross-reacts with mouse and rat. Suitable for neuroscience and axon regeneration research.

[Learn More](#)

Parameter	Value
Product Name	NgR3 Rabbit mAb
Remarks/Alias	NgR3; NGRH2
Species	Human
Gene ID (Human)	146760
Gene ID	146760
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 342-441 of human NgR3 (Q86UN2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86UN2
Protein Weight	49kDa
Shipping	Ice bag