


Anti-Nhej1 Polyclonal Antibody For Wb And Elisa - Q9H9Q4

Item Number: CM0021931



Introduction

Rabbit polyclonal antibody targeting human NHEJ1, essential for DNA double-strand break repair via NHEJ. Validated for WB and ELISA, cross-reacts with mouse. Ideal for DNA repair research.

[Learn More](#)

Parameter	Value
Product Name	NHEJ1 Rabbit pAb
Remarks/Alias	XLFI; NHEJ1
Species	Human
Gene ID (Human)	79840
Gene ID	79840
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-299 of human NHEJ1 (NP_079058.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H9Q4
Protein Weight	33kDa
Shipping	Ice bag