

Anti-Nhlh1 Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q02575

Item Number: CM0006459

Introduction

Rabbit monoclonal antibody against human NHLH1 (HEN1, NSCL, bHLHa35). Suitable for Western blot, immunofluorescence, immunoprecipitation, and ELISA. Relevant for studies of cell-type determination and neural development.

[Learn More](#)

Parameter	Value
Product Name	NHLH1 Rabbit mAb
Remarks/Alias	HEN1; NSCL; NSCL1; bHLHa35; NHLH1
Species	Human
GeneID (Human)	4807
GeneID	4807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human NHLH1 (NP_005589.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q02575
Protein Weight	15kDa
Shipping	Ice bag