

Anti-Nhlh1 Monoclonal Antibody For Wb, Ip, Elisa, Chip - Q02575

Item Number: CM0006164

Introduction

Anti-NHLH1 Rabbit monoclonal antibody validated for WB, IP, ELISA, ChIP. Targets human NHLH1 (bHLHa35), a neural basic helix-loop-helix transcription factor implicated in cell differentiation. Ideal for developmental neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	NHLH1 Rabbit mAb
Remarks/Alias	HEN1; NSCL; NSCL1; bHLHa35; NHLH1
Species	Human
GeneID (Human)	4807
GeneID	4807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human NHLH1 (NP_005589.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, ELISA, ChIP
Cross Reactivity	Human
SWISS	Q02575
Protein Weight	15kDa
Shipping	Ice bag