

Anti-Nhlh1 Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - Q02575

Item Number: CM0024139

Introduction

Rabbit polyclonal antibody targeting human NHLH1 (helix-loop-helix protein 1, bHLHa35). Validated for WB, IF/ICC, IP, and ELISA applications. Specific for human samples; ideal for developmental neurobiology research.

[Learn More](#)

Parameter	Value
Product Name	NHLH1 Rabbit pAb
Remarks/Alias	HEN1; NSCL; NSCL1; bHLHa35; NHLH1
Species	Human
GeneID (Human)	4807
GeneID	4807
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-60 of human NHLH1 (NP_005589.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human
SWISS	Q02575
Protein Weight	15kDa
Shipping	Ice bag