

Anti-Nhlrc1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q6Vvb1

Item Number: CM0027417

Introduction

Anti-NHLRC1 polyclonal antibody validated for Western blot, IHC-P, and ELISA. Reacts with human, mouse, rat. Targets E3 ubiquitin ligase malin, involved in glycogen metabolism and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	NHLRC1 Rabbit pAb
Remarks/Alias	EPM2A; EPM2B; MALIN; bA204B7.2; NHLRC1
Species	Human
GeneID (Human)	378884
GeneID	378884
Immunogen	Recombinant protein corresponding to amino acids 1-280 of human NHLRC1 (NP_940988.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6VVB1
Protein Weight	42kDa
Shipping	Ice bag