

Anti-Nhlrc1/Malin Monoclonal Antibody For Wb, If/Icc, Elisa - Q6Vvb1

Item Number: CM0008825

Introduction

Rabbit monoclonal antibody against human NHLRC1/Malin protein, validated for WB, IF/ICC, and ELISA applications. Recombinant immunogen, 42kDa predicted weight. Ideal for neurodegeneration and epilepsy research.

[Learn More](#)

Parameter	Value
Product Name	NHLRC1/Malin Rabbit mAb
Remarks/Alias	EPM2A; EPM2B; MALIN; bA204B7.2
Species	Human
Gene ID (Human)	378884
Gene ID	378884
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 79-280 of human NHLRC1 (NP_940988.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q6VVB1
Protein Weight	42kDa
Shipping	Ice bag