


Anti-Nhp2 Polyclonal Antibody For Wb, Elisa - Q9Nx24

Item Number: CM0030241



Introduction

NHP2 Rabbit Polyclonal Antibody for WB and ELISA applications. Detects human and rat NHP2 (H/ACA ribonucleoprotein complex subunit 2), involved in ribosome biogenesis and telomere maintenance. Validated in human and rat samples.

[Learn More](#)

Parameter	Value
Product Name	NHP2 Rabbit pAb
Remarks/Alias	DKCB2; NHP2P; NOLA2; NHP2
Species	Human
GeneID (Human)	55651
GeneID	55651
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human NHP2 (NP_060308.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NX24
Protein Weight	17kDa
Shipping	Ice bag