

Anti-Nid1 Rabbit Polyclonal Antibody For Wb, Elisa - P14543

Item Number: CM0030254

Introduction

NID1 Rabbit Polyclonal Antibody (pAb) validated for Western blot (WB) and ELISA in human samples. Targets nidogen-1/entactin, essential for basement membrane integrity and extracellular matrix research.

[Learn More](#)

Parameter	Value
Product Name	NID1 Rabbit pAb
Remarks/Alias	NID; NID1
Species	Human
Gene ID (Human)	4811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-400 of human NID1 (NP_002499.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P14543
Protein Weight	136kDa
Shipping	Ice bag