


Anti-Niemann Pick C1 Monoclonal Antibody For Wb - O15118

Item Number: CM0006059



Introduction

Rabbit monoclonal antibody against human Niemann Pick C1 protein (NPC1). Suitable for Western blot, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Validated for cholesterol transport research.

[Learn More](#)

Parameter	Value
Product Name	Niemann Pick C1 Rabbit mAb
Remarks/Alias	NPC; POGZ; SLC65A1; Niemann Pick C1
Species	Human
GeneID (Human)	4864
GeneID	4864
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1179-1278 of human Niemann Pick C1 (O15118).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15118
Protein Weight	142 kDa
Shipping	Ice bag