


Anti-Nifk Rabbit Polyclonal Antibody For Wb And Elisa - Q9Byg3

Item Number: CM0016773



Introduction

Rabbit polyclonal antibody targeting human NIFK (Nopp34/MKI67IP), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for nucleolar protein and cell proliferation research.

[Learn More](#)

Parameter	Value
Product Name	NIFK Rabbit pAb
Remarks/Alias	Nopp34, MKI67IP, NIFK
Species	Human
GeneID (Human)	84365
GeneID	84365
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 180-260 of human NIFK (NP_115766.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BYG3
Protein Weight	34kDa
Shipping	Ice bag