

Anti-Ninein (Nin/Sckl7) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N4C6

Item Number: CM0028844

Introduction

Rabbit polyclonal antibody against human Ninein (NIN/SCKL7). Validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for centrosome and microtubule research. SWISS: Q8N4C6.

[Learn More](#)

Parameter	Value
Product Name	NIN Rabbit pAb
Remarks/Alias	SCKL7; NIN
Species	Human
GeneID	51199
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-287 of human NIN (NP_891991.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N4C6
Protein Weight	243kDa
Shipping	Ice bag