


Anti-Ninj1 Rabbit Polyclonal Antibody For Wb And Elisa - O70131

Item Number: CM0030495



Introduction

Rabbit polyclonal antibody against mouse NINJ1 (Ninjurin-1), suitable for Western blot (WB) and ELISA. Recognizes N-terminal region (aa 1-152). Ideal for studying cell death, inflammation, and neurobiology in mouse models.

[Learn More](#)

Product Name	NINJ1 Rabbit pAb
Remarks/Alias	NIN1; NINJURIN
Species	Mouse
GeneID (human)	4814
GeneID	18081
Immunogen	Synthetic peptide (A synthetic peptide corresponding to a sequence within amino acids 1-152 of mouse NINJ1 (NP_038638.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O70131
Protein Weight	16kDa
Shipping	Ice bag