

Anti-Ninl Rabbit Polyclonal Antibody For Wb And Elisa - Q9Y2I6

Item Number: CM0028044

Introduction

Rabbit polyclonal antibody against human ninein-like protein (NINL), suitable for Western blot and ELISA applications, with cross-reactivity to mouse. Ideal for cytoskeleton and centrosome research. - Q9Y2I6

[Learn More](#)

Parameter	Value
Product Name	NINL Rabbit pAb
Remarks/Alias	NLP; NINL
Species	Human
Gene ID (Human)	22981
Gene ID	22981
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human NINL (NP_079452.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9Y2I6
Protein Weight	156kDa
Shipping	Ice bag