

Anti-Nip7 Polyclonal Antibody For Wb And Elisa - Q9Y221

Item Number: CM0031074



Introduction

NIP7 rabbit polyclonal antibody (20kDa) validated for WB and ELISA, reactive with human and mouse. Targets 60S ribosome biogenesis protein NIP7 homolog involved in nucleolar pre-rRNA processing.

[Learn More](#)

Parameter	Value
Product Name	NIP7 Rabbit pAb
Remarks/Alias	KD93; CGI-37; HSPC031; NIP7
Species	Human
Gene ID (Human)	51388
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human NIP7 (NP_057185.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y221
Protein Weight	20kDa
Shipping	Ice bag