

Anti-Nipsnap1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bpw8

Item Number: CM0019458

Introduction

Rabbit polyclonal antibody against human NIPSNAP1, validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for studying mitophagy and mitochondrial dysfunction.

[Learn More](#)

Parameter	Value
Product Name	NIPSNAP1 Rabbit pAb
Remarks/Alias	NIPSNAP1
Species	Human
GenelD (Human)	8508
GenelD	8508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 70-145 of human NIPSNAP1 (NP_003625.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BPW8
Protein Weight	33kDa
Shipping	Ice bag