

Anti-Nkd1 Rabbit Polyclonal Antibody For Wb, Elisa - Q969G9

Item Number: CM0024187

Introduction

NKD1 Rabbit Polyclonal Antibody for WB and ELISA applications, targeting human NKD1 protein. Suitable for Wnt signaling pathway research. Supplied by CamelBio for IVD raw material needs.

[Learn More](#)

Parameter	Value
Product Name	NKD1 Rabbit pAb
Remarks/Alias	Naked1; NKD1
Species	Human
GeneID (Human)	85407
GeneID	85407
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q969G9
Protein Weight	52kDa
Shipping	Ice bag