

Anti-Nkx2-1 Monoclonal Antibody For Wb And Elisa - P43699

Item Number: CM0005316

Introduction

High-quality anti-NKX2-1 rabbit monoclonal antibody validated for Western blot and ELISA. Suitable for human, mouse, rat samples. Ideal for studying thyroid transcription factor-1 (TTF-1) in thyroid and lung development.

[Learn More](#)

Parameter	Value
Product Name	NKX2-1 Rabbit mAb
Remarks/Alias	BCH; BHC; NK-2; TEBP; TTF1; NKX2A; NMTC1; T/EBP; TITF1; TTF-1; NKX2.1; NKX2-1
Species	Human
GeneID (human)	7080
GeneID	7080
Immunogen	Recombinant fusion protein containing amino acids 20-110 of human NKX2-1 (NP_003308.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43699
Protein Weight	39kDa
Shipping	Ice bag