

Anti-Nkx2-1 Rabbit Monoclonal Antibody For Ihc-P, Elisa - P43699

Item Number: CM0006658

Introduction

High-specificity Anti-NKX2-1 rabbit monoclonal antibody validated for IHC-P and ELISA, suitable for thyroid and lung research; cross-reacts with human, mouse, rat. Ideal IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	NKX2-1 Rabbit mAb
Remarks/Alias	BCH; BHC; NK-2; TEBP; TTF1; NKX2A; NMTC1; T/EBP; TITF1; TTF-1; NKX2.1
Species	Human
Gene ID (Human)	7080
Gene ID	7080
Immunogen	Recombinant protein of human NKX2-1
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43699
Protein Weight	39kDa
Shipping	Ice bag