

Anti-Nkx2-1 Rabbit Polyclonal Antibody For Wb, Elisa - P43699

Item Number: CM0030339

Introduction

CamelBio offers rabbit polyclonal antibody against NKX2-1 (TTF-1, P43699), validated for WB and ELISA. Recognizes human and rat NKX2-1, ideal for thyroid and lung development research.

[Learn More](#)

Parameter	Value
Product Name	NKX2-1 Rabbit pAb
Remarks/Alias	BCH; BHC; NK-2; TEBP; TTF1; NKX2A; NMTC1; T/EBP; TITF1; TTF-1; NKX2.1; NKX2-1
Species	Human
GeneID (Human)	7080
GeneID	7080
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human NKX2-1 (NP_003308.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P43699
Protein Weight	39kDa
Shipping	Ice bag