

Anti-Nkx2-2 Polyclonal Antibody For Wb And Elisa - 095096

Item Number: CM0030595

Introduction

NKX2-2 Rabbit Polyclonal Antibody validated for Western blot and ELISA. Reacts with human. Target is a transcriptional activator involved in pancreatic beta cell development and neuronal guidance.

[Learn More](#)

Parameter	Value
Product Name	NKX2-2 Rabbit pAb
Remarks/Alias	NKX2B; NKX2.2; NKX2-2
Species	Human
GeneID (Human)	4821
GeneID	4821
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 102-201 of human NKX2-2 (NP_002500.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O95096
Protein Weight	30kDa
Shipping	Ice bag