

Anti-Nkx2-5 Rabbit Polyclonal Antibody For Wb, Elisa - P52952

Item Number: CM0027513



Introduction

Rabbit polyclonal antibody targeting human NKX2-5 (CSX) for WB and ELISA. Cross-reacts with mouse and rat. Suitable for cardiac and spleen development research and congenital heart disease studies.

[Learn More](#)

Parameter	Value
Product Name	NKX2-5 Rabbit pAb
Remarks/Alias	CSX; CSX1; VSD3; CHNG5; HLHS2; NKX2E; NKX2.5; NKX4-1; NKX2-5
Species	Human
GeneID (Human)	1482
GeneID	1482
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-135 of human NKX2-5 (NP_004378.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52952
Protein Weight	35kDa
Shipping	Ice bag