

Anti-Nkx3-2 Rabbit Polyclonal Antibody For Wb, Elisa - P78367

Item Number: CM0029732

Introduction

Rabbit polyclonal antibody targeting human NKX3-2, suitable for WB and ELISA. Cross-reacts with mouse and rat NKX3-2. Ideal for studying chondrocyte maturation, skeletal development, and gut morphogenesis.

[Learn More](#)

Parameter	Value
Product Name	NKX3-2 Rabbit pAb
Remarks/Alias	SMMD; BAPX1; NKX3B; NKX3.2; NKX3-2
Species	Human
GeneID (Human)	579
GeneID	579
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human NKX3-2 (NP_001180.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78367
Protein Weight	35kDa
Shipping	Ice bag