

Anti-Nkx3.1 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q99801

Item Number: CM0006465



Introduction

Rabbit monoclonal antibody against NKX3.1, validated for WB, IHC-P, and ELISA. Recognizes human, mouse, and rat NKX3.1. Predicted molecular weight 26 kDa. Suitable for prostate cancer and development studies.

[Learn More](#)

Parameter	Value
Product Name	NKX3.1 Rabbit mAb
Remarks/Alias	NKX3; BAPX2; NKX3A; NKX3.1
Species	Human
GeneID	4824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human NKX3.1. (NP_006158.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99801
Protein Weight	26kDa
Shipping	Ice bag