

Anti-Nkx3.1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q99801

Item Number: CM0021049



Introduction

Rabbit polyclonal antibody against NKX3.1, validated for Western blot, immunofluorescence, and ELISA. Suitable for human and mouse samples, ideal for prostate development and cancer research.

[Learn More](#)

Parameter	Value
Product Name	NKX3.1 Rabbit pAb
Remarks/Alias	NKX3; BAPX2; NKX3A; NKX3.1
Species	Human
GeneID (Human)	4824
GeneID	4824
Immunogen	Recombinant fusion protein containing a sequence within amino acids 1-80 of human NKX3.1 (NP_006158.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99801
Protein Weight	26kDa
Shipping	Ice bag