

Anti-Nlgn4Y Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8Nfz3

Item Number: CM0029884

Introduction

Rabbit polyclonal antibody targeting NLGN4Y, validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Useful for studying cell adhesion and synapse formation in brain and testis.

[Learn More](#)

Parameter	Value
Product Name	NLGN4Y Rabbit pAb
Remarks/Alias	HNL4Y; NLGN4Y
Species	Human
Gene ID (Human)	22829
Gene ID	22829
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 122-256 of human NLGN4Y (NP_001157710.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NFZ3
Protein Weight	92kDa
Shipping	Ice bag