

Anti-NlN Polyclonal Antibody (Rabbit) For Wb And Elisa - Q9Byt8

Item Number: CM0014004



Introduction

Rabbit polyclonal antibody against human NLN (neurolysin), a mitochondrial metalloprotease. Validated for WB and ELISA; cross-reacts with human and rat. Useful for studying neuropeptide degradation and cannabinoid signaling.

[Learn More](#)

Parameter	Value
Product Name	NLN Rabbit pAb
Remarks/Alias	MEP; MOP; AGTBP; EP24.16; NLN
Species	Human
GeneID (Human)	57486
GeneID	57486
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 630-704 of human NLN (NP_065777.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BYT8
Protein Weight	81kDa
Shipping	Ice bag