


Anti-Nlrc4 Polyclonal Antibody For Wb - Q3Up24

Item Number: CM0023086



Introduction

NLRC4 rabbit polyclonal antibody for WB, IF/ICC, IF-P, ELISA. Reacts with human, mouse, rat (SWISS: Q3UP24). Suitable for inflammasome and innate immunity research.

[Learn More](#)

Parameter	Value
Product Name	NLRC4 Rabbit pAb
Remarks/Alias	CLAN; IPAF; AIFEC; CLAN1; CLANA; CLANB; CLANC; CLAND; FCAS4; CARD12; CLR2.1; NLRC4
Species	Mouse
GeneID (Human)	58484
GeneID	268973
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 71-170 of mouse NLRC4. (NP_001028539.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q3UP24
Protein Weight	116kDa
Shipping	Ice bag