


Anti-Nlrc4 Rabbit Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9Npp4

Item Number: CM0030849



Introduction

NLRC4 Rabbit Polyclonal Antibody for WB, IF/ICC, IF-P, ELISA. Validated for Human, Mouse, Rat. NLR family CARD domain-containing protein 4, key inflammasome component in innate immunity. Ideal for immunodetection and cytokine research.

[Learn More](#)

Parameter	Value
Product Name	NLRC4 Rabbit pAb
Remarks/Alias	CLAN; IPAF; AIFEC; CLAN1; CLANA; CLANB; CLANC; CLAND; FCAS4; CARD12; CLR2.1; NLRC4
Species	Human
GeneID (human)	58484
GeneID	58484
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 85-170 of human NLRC4 (NP_001186067.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NPP4
Protein Weight	116kDa
Shipping	Ice bag