

Anti-Nlrc5 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Wi3

Item Number: CM0020682

Introduction

Rabbit polyclonal antibody against human NLRC5 (NOD4, NOD27, CLR16.1), suitable for WB and ELISA with mouse cross-reactivity. Ideal for studying NF-kB/IFN signaling in innate immunity.

[Learn More](#)

Parameter	Value
Product Name	NLRC5 Rabbit pAb
Remarks/Alias	NOD4; NOD27; CLR16.1; NLRC5
Species	Human
Gene ID	84166
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human NLRC5 (NP_115582.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q86WI3
Protein Weight	205kDa
Shipping	Ice bag