


Anti-Nlrp10 Rabbit Polyclonal Antibody - Q86W26

Item Number: CM0023575



Introduction

Rabbit polyclonal antibody targeting human NLRP10 (NALP10). Suitable for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studies on innate immunity and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	NLRP10 Rabbit pAb
Remarks/Alias	NOD8, PAN5, PYNOD, NALP10, CLR11.1, NLRP10
Species	Human
Gene ID	338322
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human NLRP10 (NP_789791.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86W26
Protein Weight	75kDa
Shipping	Ice bag